

Proposal of a quality assurance program for PET-CT system



UNIVERSITY
OF TRIESTE



By Rosa A. Petit S.



I. Introduction

II. Quality assurance in PET/CT

- » QMS components
- » Professional involve
- » Importants aspects QA
- » QA and QC cycle

III. Acceptance testing

IV. List of QA/QC test


V. Other important test

VI. QC schedule

I. Introduction

3


- IAEA HUMAN HEALTH SERIES No.1- Quality assurance for PET AND PET/CT
- NEMA NU2-2007-NEMA NU2-2012
- AAPM task group No.126 – PET/CT Acceptance testing and quality assurance



II. Quality assurance in PET/CT

4

» QMS – Components



» Professionals Involve

5



» Important aspects in QA

6

Quality assurance may encompass various aspects such as quality of medical care based on specific indicators as:


- **The infection rate in the hospital.**
- **The satisfaction of patients with this care.**
- **The credential of medical staff.**
- **Any continuing education of the personal staff.**



» QA and QC cycle

7

The quality assurance and quality control cycle for a medical imaging device



III. Acceptance testing


8



Establishes the baseline performance of the equipment to which future quality tests will be compared;



Ensures that equipment (both hardware and software) performs to the manufacturer's specifications prior to final payment for the equipment;



Provides data that can give guidance in the determination of the optimal operating parameters for routine use



Ensures that the imaging equipment meets regulatory requirements for radiation safety.



NEMA NU2:


- » Spatial resolution
- » Sensitivity
- » Scatter fraction/ count rate performance
- » Image quality

ADDED TG No.126:


- » ACR (replaced image quality)
- » Uniformity
- » Alignment

IV. List of QA/QC test

10




Spatial
resolution




Sensitivity




Scatter fraction



Accuracy of
count losses



Uniformity




Alignment


V. Other important QA/QC test

11


- » Linearity
- » Low contrasts resolution
- » Artifact analysis



Linearity



Low contrast
resolution



Artifact

VI. QC schedule

12

Daily

Check singles, coincidences, timing, energy

Sinograms, tube warmup, daily air cals and daily phantom scan

Weekly

Update gains

Quarterly


Normalization and counter calibration

Annually

ACR or NEMA test, TG No.126



Pre-calibration phantom



Sample sinograms



ACR phantom images

Thanks so
much!

